

MOD.110 - ABBE Refractometer

Item Code: 43000133



DATA SHEET

Features

Brix accuracy	±0,25 %
Refractive index accuracy	±0.0002 nD

Bench-top refractometer with thermostat prism.

Description Complete with digital thermometer and calibration prism

Dimension	100 x 200 x 240 mm
Magnification of reading system	30X
Magnification of telescopic system	2X
Name	110 Abbe
Weight	2.6 kg
Brix measuring range	0,00 ... 95,00 % Brix
Refraction index measuring range	1.300 ... 1.700 nD
Temperature measuring range	0 ... 70 °C